

PRODUCT USE GUIDE FOR CLEENOL PERFUMED BACTERICIDAL MULTIPURPOSE CLEANER

Product Description

A high quality performance product providing a high standard of cleaning, degreasing and disinfection. Designed for a wide variety of applications, e.g. the cleaning of floors, walls and equipment, etc. Also for cleaning the interiors of food delivery vehicles. Its built-in residual antistatic properties protect from risk of static shock. Can be used through pressure washer or foaming equipment. Complies with BS6424 Group QAP50. Approved and tested against the MRSA bacteria.

Product Features

- An effective disinfectant, cleanser and deodoriser which cleans, disinfects and deodorises, all in one operation.
- Designed for routine use in hospitals, nursing homes, clinics, surgeries, canteens, hotels, schools and offices; for cleaning and disinfecting floors, walls, paintwork, equipment, bedframes, wipeable mattresses, sinks, drains, washrooms, lavatories and waste bins, etc.
- Complies with the European Suspension Test (BSEN1276) when diluted 1:50 (2%) after a contact time of 5 minutes, and has been found to be effective against MRSA, listeria monocytogenes, salmonella typhimurium, pseudomonas aeruginosa, escherichia coli, staphylococcus aureus and enterococcus hirae.

Directions for Use

Note: Use biocides safely. Always read the label and product information before use.

Dilution - Add 1 pump measure to 750ml pre-labelled flask.

For use in:

Washrooms: hand basins, lavatories, floors, walls

Corridors: floors, walls, skirting boards, lights.

Bedrooms: floors, walls, shelving, cabinets, commodes

Shelf Life

- It is recommended that this product is used within 12 months of purchase in the UK climate and within 6 months in warmer climates. Do not allow to freeze.

Environmental Hazards

- Neat product will cause damage to plant and animal life.
- Do not allow neat product to enter surface waterways as serious ecological damage may be caused.
- In the event of a large spillage, contact the local water authority or Environment Agency immediately.

Packaging



HD refill flask 2 x 5L

Code Nos.

083469/R 2x5L

056721/6 x6 refill flasks

Data Sheet No A19

Note: This product is only to be used for the purpose for which it was intended and in accordance with the specified instructions. This product should never be mixed with other chemical products and should be handled with care and stored carefully out of reach of children.